


## 1-Palmitoyl-2-oleoylphosphatidylethanolamine

## Chemical Properties

CAS No. :	10015-88-0
Formula:	C39H76NO8P
Molecular Weight:	718
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	1-Palmitoyl-2-oleoylphosphatidylethanolamine (POPE) is a synthetic phospholipid belonging to the PE class. It can intercalate with other amphiphilic molecules to form stable bilayers at physiological pH. 1-Palmitoyl-2-oleoylphosphatidylethanolamine is used in drug delivery research.
Targets(IC50)	Liposome

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3928 mL	6.9638 mL	13.9276 mL
5 mM	0.2786 mL	1.3928 mL	2.7855 mL
10 mM	0.1393 mL	0.6964 mL	1.3928 mL
50 mM	0.0279 mL	0.1393 mL	0.2786 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481