

## N-12:0-1-Deoxysphinganine

## Chemical Properties


CAS No. : 1246298-40-7

Formula: C<sub>30</sub>H<sub>61</sub>NO<sub>2</sub>

Molecular Weight: 467.81

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	N-12:0-1-Deoxysphinganine is a reagent used in biochemical reactions.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1376 mL	10.6881 mL	21.3762 mL
5 mM	0.4275 mL	2.1376 mL	4.2752 mL
10 mM	0.2138 mL	1.0688 mL	2.1376 mL
50 mM	0.0428 mL	0.2138 mL	0.4275 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481