


1-Palmitoyl-2-(6,7-dibromo)stearoyl-sn-glycero-3-phosphocholine

Chemical Properties

CAS No. : 324054-50-4
 Formula: C42H82Br2NO8P
 Molecular Weight: 919.88
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	1-Palmitoyl-2-(6,7-dibromo)stearoyl-sn-glycero-3-phosphocholine is a phospholipid known for its excellent biocompatibility. It can be used in the preparation of large unilamellar vesicles.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0871 mL	5.4355 mL	10.871 mL
5 mM	0.2174 mL	1.0871 mL	2.1742 mL
10 mM	0.1087 mL	0.5435 mL	1.0871 mL
50 mM	0.0217 mL	0.1087 mL	0.2174 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481