


18:1 (Δ^9 -Trans) PE

Chemical Properties

CAS No. : 19805-18-6
Formula: C₄₁H₇₈NO₈P
Molecular Weight: 744.03
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	18:1 (Δ^9 -Trans) PE is a reagent used in biochemical reactions.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.344 mL	6.7202 mL	13.4403 mL
5 mM	0.2688 mL	1.344 mL	2.6881 mL
10 mM	0.1344 mL	0.672 mL	1.344 mL
50 mM	0.0269 mL	0.1344 mL	0.2688 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481