


Phenyl 6-deoxy-2-O-(phenylmethyl)-1-thio- α -L-mannopyranoside

Chemical Properties

CAS No. : 849938-16-5
 Formula: C₁₉H₂₂O₄S
 Molecular Weight: 346.44
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Phenyl 6-deoxy-2-O-(phenylmethyl)-1-thio- α -L-mannopyranoside is a reagent used in biochemical reactions.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8865 mL	14.4325 mL	28.865 mL
5 mM	0.5773 mL	2.8865 mL	5.773 mL
10 mM	0.2887 mL	1.4433 mL	2.8865 mL
50 mM	0.0577 mL	0.2887 mL	0.5773 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481