


## 1-Palmitoyl-3-O-benzyl-rac-glycerol

## Chemical Properties

CAS No. :	1487-51-0
Formula:	C <sub>26</sub> H <sub>44</sub> O <sub>4</sub>
Molecular Weight:	420.63
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	1-Palmitoyl-3-O-benzyl-rac-glycerol is a diacylglycerol composed of Palmitic acid and Oleic acid.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3774 mL	11.8869 mL	23.7739 mL
5 mM	0.4755 mL	2.3774 mL	4.7548 mL
10 mM	0.2377 mL	1.1887 mL	2.3774 mL
50 mM	0.0475 mL	0.2377 mL	0.4755 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481