


## 2-Eicosenoic acid

## Chemical Properties

CAS No. :	109264-21-3
Formula:	C <sub>20</sub> H <sub>38</sub> O <sub>2</sub>
Molecular Weight:	310.51
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	2-Eicosenoic acid is a reagent used in biochemical reactions.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2205 mL	16.1025 mL	32.2051 mL
5 mM	0.6441 mL	3.2205 mL	6.441 mL
10 mM	0.3221 mL	1.6103 mL	3.2205 mL
50 mM	0.0644 mL	0.3221 mL	0.6441 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481