


## Diethylstilbestrol monoglucuronide

## Chemical Properties

CAS No. :	2408-40-4
Formula:	C <sub>24</sub> H <sub>28</sub> O <sub>8</sub>
Molecular Weight:	444.47
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Diethylstilbestrol monoglucuronide is a reagent used in biochemical reactions.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2499 mL	11.2494 mL	22.4987 mL
5 mM	0.450 mL	2.2499 mL	4.4997 mL
10 mM	0.225 mL	1.1249 mL	2.2499 mL
50 mM	0.045 mL	0.225 mL	0.450 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481