


## 1-Hexanoyl-2-hydroxy-sn-glycero-3-phosphocholine

## Chemical Properties

CAS No. :	58445-96-8
Formula:	C <sub>14</sub> H <sub>30</sub> NO <sub>7</sub> P
Molecular Weight:	355.36
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	1-Hexanoyl-2-hydroxy-sn-glycero-3-phosphocholine is a reagent used in biochemical reactions.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.814 mL	14.0702 mL	28.1405 mL
5 mM	0.5628 mL	2.814 mL	5.6281 mL
10 mM	0.2814 mL	1.407 mL	2.814 mL
50 mM	0.0563 mL	0.2814 mL	0.5628 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481