


D-myo-Inositol 1,3,4-tris-phosphate

Chemical Properties

CAS No. :	93133-76-7
Formula:	C ₆ H ₁₅ O ₁₅ P ₃
Molecular Weight:	420.1
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	D-myo-Inositol 1,3,4-tris-phosphate is a reagent used in biochemical reactions.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3804 mL	11.9019 mL	23.8039 mL
5 mM	0.4761 mL	2.3804 mL	4.7608 mL
10 mM	0.238 mL	1.1902 mL	2.3804 mL
50 mM	0.0476 mL	0.238 mL	0.4761 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481