


3-(Bromomethyl)phenyl isothiocyanate

Chemical Properties

CAS No. : 155863-31-3
Formula: C₈H₆BrNS
Molecular Weight: 228.11
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	3-(Bromomethyl)phenyl isothiocyanate is a reagent used in biochemical reactions.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.3838 mL	21.9192 mL	43.8385 mL
5 mM	0.8768 mL	4.3838 mL	8.7677 mL
10 mM	0.4384 mL	2.1919 mL	4.3838 mL
50 mM	0.0877 mL	0.4384 mL	0.8768 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481