


L-Idaric acid

Chemical Properties

CAS No. : 80876-58-0
 Formula: C₆H₁₀O₈
 Molecular Weight: 210.14
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	L-Idaric acid is a reagent used in biochemical reactions.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.7587 mL	23.7937 mL	47.5873 mL
5 mM	0.9517 mL	4.7587 mL	9.5175 mL
10 mM	0.4759 mL	2.3794 mL	4.7587 mL
50 mM	0.0952 mL	0.4759 mL	0.9517 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481