


2-O-(p-Nitrophenyl)- α -D-N-acetylneuraminic acid

Chemical Properties

CAS No. :	26112-88-9
Formula:	C ₁₇ H ₂₂ N ₂ O ₁₁
Molecular Weight:	430.36
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	2-O-(p-Nitrophenyl)- α -D-N-acetylneuraminic acid (PNP- α -NeuNAc) acts as a sialic acid donor in enzyme-catalyzed sialic acid transfers.
Targets(IC50)	Endogenous Metabolite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3236 mL	11.6182 mL	23.2364 mL
5 mM	0.4647 mL	2.3236 mL	4.6473 mL
10 mM	0.2324 mL	1.1618 mL	2.3236 mL
50 mM	0.0465 mL	0.2324 mL	0.4647 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481