


CBTF sodium

Chemical Properties

CAS No. : 1644038-87-8
 Formula: C₁₆H₄F₄NNaO₅S
 Molecular Weight: 421.26
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	CBTF sodium is a reagent used in biochemical reactions.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3738 mL	11.8692 mL	23.7383 mL
5 mM	0.4748 mL	2.3738 mL	4.7477 mL
10 mM	0.2374 mL	1.1869 mL	2.3738 mL
50 mM	0.0475 mL	0.2374 mL	0.4748 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481