


## 3-Aminoisonicotinic acid

## Chemical Properties

CAS No. :	7579-20-6
Formula:	C <sub>6</sub> H <sub>6</sub> N <sub>2</sub> O <sub>2</sub>
Molecular Weight:	138.12
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	3-Aminoisonicotinic acid is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.2401 mL	36.2004 mL	72.4008 mL
5 mM	1.448 mL	7.2401 mL	14.4802 mL
10 mM	0.724 mL	3.620 mL	7.2401 mL
50 mM	0.1448 mL	0.724 mL	1.448 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481