


Nonenylsuccinic anhydride

Chemical Properties

CAS No. : 28928-97-4
 Formula: C₁₃H₂₀O₃
 Molecular Weight: 224.3
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Nonenylsuccinic anhydride is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.4583 mL	22.2916 mL	44.5831 mL
5 mM	0.8917 mL	4.4583 mL	8.9166 mL
10 mM	0.4458 mL	2.2292 mL	4.4583 mL
50 mM	0.0892 mL	0.4458 mL	0.8917 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481