


Pentamethoxy red

Chemical Properties

CAS No. : 1755-51-7
 Formula: C₂₄H₂₆O₆
 Molecular Weight: 410.46
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Pentamethoxy red is a reagent used in biochemical reactions.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4363 mL	12.1815 mL	24.3629 mL
5 mM	0.4873 mL	2.4363 mL	4.8726 mL
10 mM	0.2436 mL	1.2181 mL	2.4363 mL
50 mM	0.0487 mL	0.2436 mL	0.4873 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481