


2,2-Diethoxyacetophenone

Chemical Properties

CAS No. :	6175-45-7
Formula:	C ₁₂ H ₁₆ O ₃
Molecular Weight:	208.25
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2,2-Diethoxyacetophenone (Photoinitiator-DEAP) is a photoinitiator useful for preparing hydrogel adsorbents.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.8019 mL	24.0096 mL	48.0192 mL
5 mM	0.9604 mL	4.8019 mL	9.6038 mL
10 mM	0.4802 mL	2.401 mL	4.8019 mL
50 mM	0.096 mL	0.4802 mL	0.9604 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481