


## 1-Naphthyl-β-D-glucuronide sodium

### Chemical Properties

CAS No. : 83833-12-9  
 Formula: C<sub>16</sub>H<sub>15</sub>NaO<sub>7</sub>  
 Molecular Weight: 342.28  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
*Actual storage temperature shall be subject to the COA.*



### Biological Description

Description	1-Naphthyl-β-D-glucuronide sodium is a reagent used in biochemical reactions.
Targets(IC50)	Others

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9216 mL	14.6079 mL	29.2158 mL
5 mM	0.5843 mL	2.9216 mL	5.8432 mL
10 mM	0.2922 mL	1.4608 mL	2.9216 mL
50 mM	0.0584 mL	0.2922 mL	0.5843 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481