


1-desoxymethylsphinganine

Chemical Properties

CAS No. :	1219484-98-6
Formula:	C17H37NO
Molecular Weight:	271.48
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	1-Desoxymethylsphinganine is a sphingoid base and the 1-deoxymethyl derivative of sphinganine.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6835 mL	18.4176 mL	36.8351 mL
5 mM	0.7367 mL	3.6835 mL	7.367 mL
10 mM	0.3684 mL	1.8418 mL	3.6835 mL
50 mM	0.0737 mL	0.3684 mL	0.7367 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481