


D-Ribopyranosylamine

Chemical Properties

CAS No. :	43179-09-5
Formula:	C ₅ H ₁₁ NO ₄
Molecular Weight:	149.15
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	D-Ribopyranosylamine is a biochemical reagent used in glyobiology research, which focuses on the study of sugar structures, synthesis, biology, and evolution. This field involves carbohydrate chemistry, glycoconjugate synthesis and degradation enzymology, protein-glycan interactions, and the role of glycans in biological systems. Glyobiology is closely connected with fundamental research, biomedicine, and biotechnology.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.7047 mL	33.5233 mL	67.0466 mL
5 mM	1.3409 mL	6.7047 mL	13.4093 mL
10 mM	0.6705 mL	3.3523 mL	6.7047 mL
50 mM	0.1341 mL	0.6705 mL	1.3409 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481