


meso-1,2-Bis(4-methoxyphenyl)ethylenediamine

Chemical Properties

CAS No. : 117106-36-2
 Formula: C₁₆H₂₀N₂O₂
 Molecular Weight: 272.34
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	meso-1,2-Bis(4-methoxyphenyl)ethylenediamine is a reagent used in biochemical reactions.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6719 mL	18.3594 mL	36.7188 mL
5 mM	0.7344 mL	3.6719 mL	7.3438 mL
10 mM	0.3672 mL	1.8359 mL	3.6719 mL
50 mM	0.0734 mL	0.3672 mL	0.7344 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481