


Methyl cis-11-octadecenoate

Chemical Properties

CAS No. : 1937-63-9
 Formula: C₁₉H₃₆O₂
 Molecular Weight: 296.49
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Methyl cis-11-octadecenoate is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3728 mL	16.864 mL	33.728 mL
5 mM	0.6746 mL	3.3728 mL	6.7456 mL
10 mM	0.3373 mL	1.6864 mL	3.3728 mL
50 mM	0.0675 mL	0.3373 mL	0.6746 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481