

Stigmasta-3,5-diene

Chemical Properties


CAS No. : 4970-37-0

Formula: C₂₉H₄₈

Molecular Weight: 396.69

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Stigmasta-3,5-diene is a reagent used in biochemical reactions.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5209 mL	12.6043 mL	25.2086 mL
5 mM	0.5042 mL	2.5209 mL	5.0417 mL
10 mM	0.2521 mL	1.2604 mL	2.5209 mL
50 mM	0.0504 mL	0.2521 mL	0.5042 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481