


## 2-Hydroxydodecanoic acid

## Chemical Properties

CAS No. :	2984-55-6
Formula:	C <sub>12</sub> H <sub>24</sub> O <sub>3</sub>
Molecular Weight:	216.32
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	2-Hydroxydodecanoic acid is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.6228 mL	23.1139 mL	46.2278 mL
5 mM	0.9246 mL	4.6228 mL	9.2456 mL
10 mM	0.4623 mL	2.3114 mL	4.6228 mL
50 mM	0.0925 mL	0.4623 mL	0.9246 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481