


## N-d-Biotinyl-7-amino-4-methylcoumarin

### Chemical Properties

CAS No. :	191223-35-5
Formula:	C <sub>20</sub> H <sub>23</sub> N <sub>3</sub> O <sub>4</sub> S
Molecular Weight:	401.48
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



### Biological Description

Description	N-d-Biotinyl-7-amino-4-methylcoumarin is a reagent used in biochemical reactions.
Targets(IC50)	Others

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4908 mL	12.4539 mL	24.9078 mL
5 mM	0.4982 mL	2.4908 mL	4.9816 mL
10 mM	0.2491 mL	1.2454 mL	2.4908 mL
50 mM	0.0498 mL	0.2491 mL	0.4982 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481