

## Rifampin-BSA

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Rifampin-BSA is a hapten-carrier protein conjugate formed by binding Rifampin to BSA (bovine serum albumin). It enhances Rifampin-induced immune responses by combining a hapten with a protein adjuvant, commonly used in vaccine or antibody preparation studies.
Targets(IC50)	Antibiotic

## Reference

Weiss AM, et al. Site-specific antigen-adjuvant conjugation using cell-free protein synthesis enhances antigen presentation and CD8+ T-cell response. Sci Rep. 2021 Mar 18;11(1):6267.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481