

## Pyrazinamide-KLH

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Pyrazinamide-KLH is a hapten-carrier protein conjugate formed by coupling Pyrazinamide to keyhole limpet hemocyanin (KLH). As a hapten, Pyrazinamide alone cannot induce an immune response and must be conjugated to a large molecular protein to become a complete antigen with immunogenicity. Pyrazinamide-KLH may be supplied as a lyophilized powder or as a solution. The solution is prepared in 0.01M PBS (pH 7.4) buffer at a typical concentration of 1mg/mL.
Targets(IC50)	Antibacterial,Antibiotic,Autophagy

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481