

## Periplogenin-KLH

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

#### Description

Periplogenin-KLH is a hapten-carrier protein conjugate formed by coupling Periplogenin to keyhole limpet hemocyanin (KLH). As a hapten, Periplogenin alone cannot induce an immune response and must be conjugated to a large molecular protein to become a complete antigen with immunogenicity. Periplogenin-KLH may be supplied as a lyophilized powder or as a solution. The solution is prepared in 0.01M TBS (pH 7.4) buffer at a typical concentration of 1mg/mL.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481