

AOZ-BSA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	AOZ-BSA is a hapten-carrier protein conjugate formed by coupling AOZ to bovine serum albumin (BSA). As a hapten, AOZ alone cannot induce an immune response and must be conjugated to a large molecular protein to become a complete antigen with immunogenicity.
Targets(IC50)	Drug Metabolite

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481