

Angiotensin II-BSA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Angiotensin II-BSA is a hapten-carrier protein conjugate formed by coupling Angiotensin II to bovine serum albumin (BSA). As a hapten, Angiotensin II alone cannot induce an immune response and must be conjugated to a large molecular protein to become a complete antigen with immunogenicity. The sequence of Angiotensin II-BSA is DRVYIHPF-BSA.
Targets(IC50)	RAAS

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481