

## Desthiobiotin-PEG3-triazole-DBCO-sulfo-Maleimide

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Desthiobiotin-PEG3-triazole-DBCO-sulfo-Maleimide serves as an ADC linker, utilized in the synthesis of antibody-drug conjugates (ADCs). This compound is a click chemistry reagent, featuring a DBCO group that enables strain-promoted alkyne-azide cycloaddition (SPAAC) with molecules containing azide groups.
Targets(IC50)	ADC Linker

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481