

Sodium lactate ringer's solution

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

4°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Sodium lactate ringer's solution is an isotonic crystalloid solution used to regulate body fluids, electrolytes, and acid-base balance, and it is effective for fluid replenishment during the perioperative period (perioperative). This solution consists of a sterile aqueous solution of sodium lactate, sodium chloride, potassium chloride, and calcium chloride (RL, 130 mEq sodium/L, 274 mOsm/L).
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481