

L-Histidine buffer (Sterile)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Store at 4°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	L-Histidine buffer (Sterile) is a widely used antibody-storage formulation containing 20 mM L-Histidine, 240 mM sucrose, 150 mM NaCl, and 0.05% Tween-80 at pH 5.8. L-Histidine buffer (Sterile) supports long-term protein stability by providing osmoprotection, isotonicity, and controlled acidity, making it suitable for monoclonal antibody preservation and biologics development.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481