

Griffonia Simplicifolia Lectin I Isolectin B4 (Biotinylated)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Griffonia Simplicifolia Lectin I Isolectin B4 (Biotinylated) is a plant-derived lectin that can specifically interact with biomolecules such as polysaccharides and peptides, and is therefore used as a probe tool in biological analysis. Griffonia Simplicifolia Lectin I Isolectin B4 (Biotinylated) is also a substance that can be used as a biomaterial or organic compound in life science research.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481