

## Bovine Serum Albumin-PE

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Bovine Serum Albumin-PE is a biochemical reagent conjugated with the fluorophore PE (Ex is 565 nm; Em is 578 nm). This compound is utilized for labeling or detecting specific antigens, molecules, or cell structures.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481