

Ovalbumin-AF647

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Ovalbumin-AF647 is an egg white protein conjugated with the far-red fluorescent dye AF647 (Ex=650 nm; Em=665 nm). It serves as a fluorescent marker for labeling or detecting specific antigens, molecules, or cellular structures.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481