

## Datura Stramonium Lectin

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

|               |  |
|---------------|--|
| Description   | Datura Stramonium Lectin is a plant-derived lectin that can specifically interact with biomolecules such as polysaccharides and peptides, and is therefore used as a probe tool in biological analysis. Datura Stramonium Lectin is also a substance that can be used as a biomaterial or organic compound in life science research. |
| Targets(IC50) | Others   |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481