

Type III collagen (from bovine skin)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Type III collagen (from bovine skin) is a collagen protein capable of stimulating the immune system in animals to generate specific antibodies against it. Additionally, this collagen serves as a hydrolytic substrate for matrix metalloproteinases (MMPs).
Targets(IC50)	MMP

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481