

## Type XI collagen (from bovine articular cartilage)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Type XI collagen (from bovine articular cartilage) is an immunization-grade collagen protein that induces the animal immune system to produce specific antibodies against it. Additionally, this collagen serves as a hydrolytic substrate for matrix metalloproteinases (MMP).
Targets(IC50)	MMP

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481