

## Type I collagen (from mouse skin)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Type I collagen (from mouse skin) is a highly purified form of type I collagen derived from mouse skin. It is capable of stimulating the immune system in animals to produce specific antibodies against this collagen. Additionally, collagen serves as a hydrolytic substrate for matrix metalloproteinases (MMPs).
Targets(IC50)	MMP

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481