

Type II collagen (from mouse sternal cartilage)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Type II collagen (from mouse sternal cartilage) is designed for immunization purposes. It can stimulate an animal's immune system to produce specific antibodies against this collagen. Additionally, collagen serves as a hydrolytic substrate for MMP.
Targets(IC50)	MMP

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481