

Type V collagen (from mouse intestine)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Type V collagen (from mouse intestine) serves as an immunization-grade protein that can elicit a specific immune response in animals, generating antibodies directed against this collagen. Additionally, this collagen functions as a hydrolytic substrate for matrix metalloproteinases (MMPs).
Targets(IC50)	MMP

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481