

Type XI collagen (from porcine articular cartilage)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Type XI collagen (from porcine articular cartilage) serves as a collagen protein that can stimulate the animal immune system to produce specific antibodies against this collagen. Additionally, the collagen protein acts as a hydrolytic substrate for MMP.
Targets(IC50)	MMP

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481