

Type I Collagen (from Mouse)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Type I Collagen, T Cell Grade, from Mouse possesses unique immunological properties and plays a significant role in maintaining the structure and function of skin, bones, and tendons. Additionally, it can serve as a stimulatory antigen for T cells in in vitro culture systems.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481