

## Methamphetamine-BSA

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Methamphetamine-BSA is a hapten-carrier protein conjugate formed by coupling Methamphetamine to bovine serum albumin (BSA). As a hapten, Methamphetamine alone cannot induce an immune response and must be conjugated to a large molecular protein to become a complete antigen with immunogenicity. Methamphetamine-BSA may be supplied as a lyophilized powder or as a solution. The solution is prepared in 0.01M PBS (pH 7.4) buffer at a typical concentration of 1mg/mL.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481