

Bismuth Sulfite Agar Medium

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|--|
| Description | Bismuth Sulfite Agar Medium is utilized for the selective isolation and cultivation of Salmonella. The composition of Bismuth Sulfite Agar Medium includes tryptone, peptone, beef extract, ferrous sulfate, bismuth sulfite, disodium phosphate, dextrose, brilliant green, and agar. |
| Targets(IC50) | Others |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481