

Blood Enrichment Medium

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Blood Enrichment Medium is utilized for the enrichment culture of pathogens present in the blood. Its composition includes peptone, beef extract, sodium chloride, glucose, sodium citrate, p-aminobenzoic acid, magnesium sulfate, and phenol red.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481