

Casein Agar

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Casein Agar is applicable for the casein decomposition test of Bacillus cereus. Each liter of Casein Agar medium contains: 10.0 g casein, 3.0 g beef extract powder, 5.0 g sodium chloride, 2.0 g dipotassium phosphate, and 15.0 g agar.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481