

ST3 β -Gal α -2,3-Sialyltransferase 2

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	ST3 beta-Gal alpha-2,3-Sialyltransferase 2 (EC:2.4.3.4, ST3GAL2) catalyzes the transfer of sialic acid from CMP-sialic acid to galactose-containing substrates.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481